## MODEL CLD-500 CLASS 100 CLEAN ROOM CABINET OVEN, Serial No. 730047

for sterilizing glassware

- -maximum operating temperature: 500F
- -work space dimensions: 36" wide x 36" deep x 39" high
- -20 KW installed in Incoloy sheathed tubular heating elements
- -1200 CFM, 1-1/2-HP recirculating blower providing horizontal airflow
- -4" thick insulated walls
- -type 304, 2B finish stainless steel interior with 1/2" inside radius corners, seams welded, ground and polished
- -stainless steel exterior with #4 brushed finish
- -door front and rear
- -wall sealing flange
- -red pilot light to indicate when opposite door is open
- -solenoid operated door locks
- -chrome plated door hardware
- -12" x 12" x 6" thick HEPA fresh air filter with 2" prefilter
- -30" x 24" x 6" thick stainless steel high temperature HEPA recirculating filter
- -minihelic pressure gauge and DOP validation port across each set of filters
- -digital programming temperature controller
- -manual reset excess temperature controller with separate contactors
- -recirculating blower air flow safety switch
- -10" diameter circular chart temperature recorder
- -SCR power controller
- -stainless steel loading truck with 6 removable pans
- -transfer dolly including oven-to-dolly latch and dolly-to-truck latch
- -fused disconnect switch

