

**MODIFIED MODEL HA-500 HORIZONTAL AIR FLOW CABINET OVEN,  
Serial No. 86942D0406 for laboratory testing**

- maximum operating temperature: 500 F
- work space dimensions: 20" wide x 20" deep x 20" high
- 8 KW installed in Incoloy sheathed tubular heating elements
- 400 CFM, 1/3-HP recirculating blower providing horizontal air flow
- 4" insulated walls
- type 304, 2B finish stainless steel interior with continuously back-welded seams
- stainless steel exterior with #4 brushed finish
- safety equipment for handling flammable solvents including explosion venting door hardware
- digital indicating temperature controller



**THE GRIEVE CORPORATION**

500 Hart Road, Round Lake, Illinois 60073-1098 USA  
(847) 546-8225 Fax: (847) 546-9210

**No.1002**