

MODIFIED MODEL LO-201C LABORATORY OVEN, Serial No. 630392

for preheating 15 lb. molds

- maximum operating temperature: 200C
- work space dimensions: 12" wide x 10" deep x 10" high
- 800 watts installed in an Incoloy sheathed tubular heating element
- gravity convection

- 1" insulated walls
- aluminized steel interior
- stainless steel exterior with #4 brushed finish
- heavy duty door latch
- reinforced door hinge
- shelf supports bolted through walls to support load at double shelf
- temperature control shaft recessed and slotted to avoid set point change



THE GRIEVE CORPORATION

500 Hart Road, Round Lake, Illinois 60073-2898 USA
(847) 546-8225 Fax: (847) 546-9210

No.969