


**MODIFIED MODEL AB-650 UNIVERSAL OVEN, Serial No. 350661**

for heat depyrogenation

- maximum operating temperature: 650°F
- workspace dimensions: 36" wide x 36" deep x 36" high
- 12 KW installed in Incoloy sheathed tubular heating elements
- 600 CFM, 1/2-HP recirculating blower providing universal airflow
- 6" insulated walls
- aluminized steel exterior
- Type 304, 2B finish stainless steel interior
- positive latching door hardware
- solenoid operated door lock
- safety equipment for handling flammable solvents (including explosion venting panel in floor of oven)
- integral leg stand



**THE GRIEVE CORPORATION**

500 Hart Road, Round Lake, Illinois 60073-2898 USA  
(847) 546-8225 Fax: (847) 546-9210

**No.898**